

METRO TECHPACK SYSTEM

32600-Series

Available in Impervious UP Fabric:GNUP (Green), RBUP (Red), ORUP (Orange), BKUP (Black) or RDUP (Red)

Built for Life. In the USA.

Compact. Protocol specific. Fast deployment, simple decon and easy restocking.

Rugged, durable and backed by a lifetime warranty on materials and workmanship.



Flexible carry options; padded tablet pocket



Huge fold-down first-in pocket with exterior O2 access



Deluxe, UP Pack Straps for comfort, injury prevention and easy decontamination



Pack straps have a two point adjustment system for a secure fit whether climbing, running, walking or jumping

Choose from BLS, Advanced Airway and ALS Setups

32600-02-AdvAir Setup



32600-ALS Setup

BLS Page



32600-02-BLS Setup



BLS MedModule

Optional O2 Module



Optional AED Case



Pack straps and removable waist strap store inside the ergonomically padded, zippered rear compartment

360° reflective

Semi-rigid, padded shell is lined with high visibility safety yellow for enhanced visibility.

No hook and loop in the Pack, Pages or Modules for easy decon inside and out

Item #	Description	Dimensions
ID: 32600-02-BLS BT: 2522-32602	Metro TechPack O2 BLS Setup includes: <ul style="list-style-type: none"> MTP O2 Case (O2 strap; carries one Page) BLS Page (Bleeding and burn control) BLS MedModule (BLS meds and delivery) 	20.5" H x 14" W x 9" D (Wt: 5.7 lbs) 14" H x 13.5" W x 4" D (Wt: 2.0 lbs)
ID: 32600-02-AdvAIR BT: 2522-32600	Metro TechPack O2 Advanced Airway Setup: <ul style="list-style-type: none"> MTP O2 Case (O2 strap; carries one Page) Adv Air Page (Adult/Pediatric intubation) 	20.5" H x 14" W x 9" D (Wt: 5.7 lbs) 14" H x 13.5" W x 4" D (Wt: 2.0 lbs)
ID: 32600-ALS BT: 2521-32601	Metro TechPack ALS Setup includes: <ul style="list-style-type: none"> MTP Case (No O2 strap; carries two Pages) MTP IV/Infusion Page (IO, IV, Infusion) MTP ALS MedModule Accepts optional 2nd page 	20.5" H x 14" W x 9" D (Wt: 5.7 lbs) 14" H x 13.5" W x 4" D (Wt: 2.0 lbs)
32610-UP	Optional MTP O2 Wrap (O2 transport w/chair or	20.5" H x 14" W x 9" D
32615-UP	Optional AED Briefcase	14" H x 13.25" W x 4.4" D



CUSTOM LOGO
AVAILABLE
see details on website

PACK ONE

WWW.IRONDUCK.COM



IRON DUCK
BUILT FOR LIFE

METRO TECHPACK SYSTEM

32600-Series



32600-02-BLS Setup and 32600-02 Advanced Airway Setup

Organize basic airway supplies in the chassis and lid pockets of the O2 Backpack:

- Pediatric sized nebulizers, airways and non-re-breathers in one long lid pocket and adult sizes in the other.
- Tank straps w/quick release buckles hold a "D" or "Jumbo D" tank, with or without the optional O2 Wrap
- Adult and pediatric sizes of BVMs and Nebulizers fit to the side of O2
- All pockets constructed of durable, non-yellowing, clear vinyl for a crisp, clear view of equipment within
- Attach the BLS Page to organize Bleeding control on one side, Burn care on the other
- Or Attach the Advanced Airway Page to configure as Advanced Air Setup; organize Pediatric sized blades, Macgill forceps, ET tubes and tube holders on the front, Adult on the back

32600-ALS Setup

Organize basic trauma supplies in the chassis and first in diagnostics in the fold-down lid pockets:

- Bleeding control in one long pocket, burn treatments in the other
- Organize IO, IV and Infusion supplies in the IV/ Infusion Page; IO and accessories in the pouch, Micro and Macro Setups in the long center pocket, syringes at lower right and needle kits lower left
- Two sets of magnetic hardware mean a second page (available at additional cost) will store Advanced Airway supplies - Adult on one side, and Pediatric on the other
- All pockets constructed of durable, non-yellowing, clear vinyl for a crisp, clear view of equipment within for years to come
- The ALS MedModule stores a full complement of meds and delivery devices in the removable, color coded pages

Huge fold-down, first-in pocket holds two sizes of stethoscopes, BP cuffs, gloves, flashlight, trauma sheers, glucometer kit, gloves, wipes and mask.



Monitor regulator levels through the same fold-down first-in pocket.

Zippered, padded exterior pocket holds any sized tablet or document.

The briefcase-style modules function like Pages but are easily removed and quickly reattached with magnetic hardware. Use dual zippers with equipment seals. Stock extra pages and swap as needed.



**PACK ONE
PACK ONE
PACK ONE**

WWW.IRONDUCK.COM

20 VETERANS DRIVE
CHICOPEE MA 01022
413-593-3300

