

FLAT BOARD



Part# 35841

Pedi Air Align Spineboard

- This pediatric sized board is especially designed for children infant to 5 years old.
- The board's patented Dual Plane Head Drop System supports the larger occipital area of a child's head and allows you to immobilize a patient, establish neutral spine alignment and maintain their airway all at the same time.
- The only pediatric immobilization board that is roto-molded from high-density polyethylene (HDPE) to create a one piece seamless product, eliminating cross contamination.
- X-Ray and CT Scan translucent, MRI compatible.

COLORS:

- Yellow • Blue • Pink • Teal

OPTIONS:

- Custom graphics available. Additional fees apply.



Part# 35840

Pedi Air Complete

Kit Contains:

- 1 Pedi Board
- 1 Set Restraint Straps (Pedi)
- 1 Shoulder Strap
- 1 Chest Strap
- 4 Body Straps
- 2 Wrist Restraints
- 6 Disposable Head Blocks



• **Size:** 48" x 12" x 1.5"

• **Weight:** 7.5 lbs.

• **Capacity:** 1lb. single patient



BUILT IN AMERICA



413-377-0030
www.ironduck.com

