




Midwife BIRTHING BAGS



Durable, versatile bags for wherever your work takes you



34016 D

Midwife Breathsaver

Dimensions: 27"L x 13"W x 11"H
3861 cu in

- Main compartment holds a D or Jumbo D oxygen tank with straps to secure.
- Compartments have mesh pouches and elastic tiedowns to organize oxygen supplies.
- Includes shoulder strap.
- Interior liner is removable for cleaning.



#32350

Midwife Prenatal Case

Dimensions: 19"L x 14"W x 12"H
3132 cu in

- Main compartment is divided into three sections with large zippered pockets and elastic tool holders.
- Front and rear pockets have multiple interior pouches and elastic holders.
- Heavy duty handles wrap around bag for even weight distribution.
- Includes shoulder strap.



32325

Midwife Ultra Sofbox Plus

Dimensions: 23"L x 17"W x 12"H
4692 cu in

- Main compartment's moveable dividers configure to your needs.
- Top lid has clear mesh reinforced pocket and elastic tool holders.
- Large pockets on each end provide storage without bulk.
- Impact resistant runners; UP material protects bottom
- Deluxe ergonomic padded shoulder strap included.
- Interior liner is removable for cleaning.



#32470

Midwife Backpack Plus

Dimensions: 23"L x 14"W x 12"D

- Front bottom pouch can hold an AED, which can be secured with quick release tie downs.
- Main compartment offers multiple customizable storage options
- Removable organizational pouch has flat zippered and clear gusseted pockets.
- O₂ compartment has quick-release tie-downs to secure a "D" or a Jumbo "D" oxygen tank.
- Corners are reinforced with high density plastic to prevent wear.
- Bag can be carried using the ergonomic backpack straps, the side, or top molded handles.

To view more product photos, go to:
www.ironduck.com/store/midwife-bags



www.ironduck.com • info@ironduck.com • 413-377-0030



BUILT IN AMERICA