

# RAPID RESPONSE PACK

32670-XX\*



## \*Available Colors

Non-standard colors available – extended lead time for production.

1,000 Denier Nylon: (XX = RD), (XX = RB), (XX = OR), (XX = GN), (XX = BK)

Universal Precautions: (XX = RDUP), (XX = RBUP), (XX = ORUP), (XX = GNUP), (XX = BKUP)

SALES: 1-413-593-3300x101  
FAX: 1-413-593-5800



FRONT

BACK



SIDE

SIDE

PROFILE



BACKPACK STRAPS

## EXTERIOR FEATURES



TOP

Iron Duck's Rapid Response Pack is versatile enough to fit the needs of First-In ALS, Rapid Response, USAR and/or Special Event Teams. Built to last, in the USA, this rugged, compact backpack is efficiently organized, quick to deploy and easy to decontaminate.

- A** Available in impervious Universal Precautions fabric or 1000D Nylon.
- B** Slim, adjustable, ergonomic backpack straps feature extra padding at the clavical and D-rings for the attachment of small modules. Backpack Straps fold up compactly beneath the padded lower back support with magnetic closures. Removable waist strap.
- C** Carry with riveted molded rubber top or side drag handles.
- D** Clear vinyl, business card sized ID holder on top and or customized screen print for easy identification.

**Built for Life. In the USA.**



**PACK ONE**

Item #	Product	Exterior Dimensions
32670-XXUP	Rapid Response Pack	18"Hx ***12"W x 10"D (4 lbs ?) Depth includes tucked straps
Pocket	Features	Exterior Dimensions
E	Fold-down, first-in pocket <ul style="list-style-type: none"> <li>Holds first-in diagnostic supplies; blood pressure cuffs, stethoscoped, gloves, glucose kits</li> <li>Dual zippers with large, hi-viz pull tabs</li> </ul>	6"H x 12"W x 2"D
F	Bottom exterior pocket <ul style="list-style-type: none"> <li>Inside lid features classic trauma set-up. Divided clear vinyl pocket features three compartments with easy access via one hood with magnetic closure on either side. A row of elastic tool holders over a flat zippered pocket completes the bottom exterior pocket lid.</li> <li>Opposite, the exterior pocket main compartment features one clear vinyl gusset; hood closes with magnetic closure</li> </ul>	12"H x 12"W x 2"D Each gusset: 6"H x 4"W x 2"D  10"H x 12"W x 2"D


# RAPID RESPONSE PACK

## 32670-XXUP



### INTERIOR FEATURES

- A** First-in, fold-down pocket holds diagnostic supplies: Blood Pressure Cuff, Glucose, Gloves and Mask
- B** The lid of the bottom exterior pocket features the classic trauma design; a flat pocket w/dual zipper closure behind a row of elastic tool holders and gusset pocket divided into three compartments
- C** Opposite, a single clear vinyl pocket with magnetic closures on hood, is large enough to hold compact AEDs
- D** The main compartment lid features
- E** **Optional** Utility Module
- F** **Optional** Divider System



Pocket	Features	Interior Dimensions
A	Fold-down, First-in Pocket <ul style="list-style-type: none"> <li>Holds gloves, mask and first-in diagnostic equipment (one size blood pressure cuff, stethoscope, glucose kit)</li> </ul>	6"H x 12"W x 2"D
B	Bottom Exterior Pocket - Inside Lid <ul style="list-style-type: none"> <li>Flat pocket w/dual zipper closure</li> <li>Row of 1" elastic tool holders</li> <li>Gusseted, clear vinyl pocket w/structured hood and magnetic closures, divided into two small and one large compartment</li> </ul>	11.5"L x 11.5"W  Two small: 6"H x 2.5"W x 2.5"D One large: 6"H x 6"W x 2.5"D
C	Bottom Exterior Pocket - Main Compartment <ul style="list-style-type: none"> <li>Clear vinyl pocket w/structured hood and magnetic closure holds most small AEDs</li> </ul>	11.5"H x 11.5"W x 2.5"D
D	Main Compartment - Inside Lid <ul style="list-style-type: none"> <li>Just under the Fold-down, First-in Pocket are two long, clear vinyl pockets each w/dual zipper closures</li> </ul>	11.5"H x 5.5"W x 2"D
E	<b>Optional Utility Module</b> <ul style="list-style-type: none"> <li>Two adjustable, dividers w/clear vinyl pockets on either side</li> <li>Fold out lid with two clear vinyl, elastic topped pockets</li> </ul>	11.5"L x 4"W x 4"D
F	<b>Optional Adjustable Divider System w/Clear Vinyl Lid and Magnetic Closure</b> <ul style="list-style-type: none"> <li>Move or remove dividers to fit almost any configuration.</li> </ul>	11.5"H x 11.5"W x 4"D