

# Ultra 2020 Spineboard

Item #35740

## A Bigger, Wider, Spineboard that FLOATS.

Perfect for pools, parks and beaches where safety is a concern

**NEW**



### Features:

- Ideal for water rescue; helps to float victims; Meets **OSHA EM385-1-1 2014 Code** for Diving Safety
- Wider board accommodates wider patients
- 14 extra-large, non slip textured handles help ensure safe handling in water
- Concave patient surface provides patient comfort and support while minimizing lateral movement
- Spread out runners: wider in torso region (crucial for added support)
- Compatible with most standard head immobilizers and restraints
- Multiple strap holes; available with/without 12 speed clip pins (*Strapping pins not recommended with water use*)
- CT MRI and X-ray compatible
- Fits most vehicle storage areas
- Manufactured in the USA to USA Quality Standards
- Colors: Yellow, Blue. *Contact your dealer for other colors*

### Specifications:

- Size: 74"L x 18.5"W x 3.5"H
- Board Weight: 14 lbs.
- Weight Capacity: 1000 lbs.
- Float Capacity: 200 lbs spread over surface of board





**413-377-0030**  
**www.ironduck.com** • info@ironduck.com  
**20 Veterans Drive • Chicopee, Massachusetts**

**Be Ready.**  
**Be Prepared.**

